

Warning:

This product is for research use only.

Not recommended or intended for diagnosis of disease in humans or animals. Do not use internally or externally in humans or animals.

Contains sodium azide. For further instructions see material safety data sheet (MSDS).

Introduction

Streptavidin is a biotin-binding protein, which is excreted by the bacterium *Streptomyces avidinii*. Streptavidin is a non-glycosylated 55 kDa protein, which is composed of 4 essentially identical polypeptide chains. Unlike avidin, with an isoelectric point of 10.5, streptavidin is neutral at pH 7. This is the reason, why streptavidin shows considerably less non-specific binding and thus less background as avidin, which makes it the reagent of choice for the detection of biotinylated biomolecules.

Formulation

Our Oyster[®]-Streptavidin conjugates are supplied as lyophilized powder in 1 mg quantities. For stabilization and preservation the conjugate is lyophilized from a solution containing 10 % bovine serum albumin and 0.05 % sodium azide in phosphate buffered saline (PBS), pH 7.4. Highly purified and electrophoretically homogeneous streptavidin was conjugated with Oyster[®] dyes under optimized conditions and is free of non-conjugated fluorophore.

Storage

Avoid freeze-thaw cycles. Protect from light during storage. Keep the container at 4 °C before opening. For long term storage, addition of 40 % glycerol (w/v) and storage at -20 °C is recommended. Under these conditions the product is stable for at least 6 months.

Applications

Fluorescent streptavidin conjugates can be used in many common applications including immune-histochemistry, Western blotting, Fluorescence *In Situ* Hybridisation (FISH) and ELISA.

The exact working concentrations should be determined individually by the researcher, as the appropriate concentration is affected by the concentration of the biotinylated biomolecule, the incubation time or the affinity of the antigens or antibodies used.

We recommend the following concentrations to start with:

Immune-histochemistry:	1:200 to 1:500
Western blot	1:1,000 to 1:2,000
ELISA	1:2,000 to 1:5,000

Oyster[®]-Streptavidin conjugates are available as

Label	Ex_{max}	Em_{max}	Product code
Oyster [®] 550	555	574	OY-550-SAV-1x1
Oyster [®] 650	655	674	OY-650-SAV-1x1

Oyster[®] 550-Streptavidin

The Oyster[®]-550-Streptavidin is a bright, orange fluorescent conjugate, which is excited best by the 546 nm line of a mercury-arc lamp or the 543 nm of a He-Ne laser. Oyster[®]-550 is spectrally similar to Cy3[™], Tetramethylrhodamine (TRITC) or Alexa Fluor[®] 555.

Oyster[®] 650-Streptavidin

The Oyster[®]-650-Streptavidin is a bright, red fluorescent conjugate, which is excited best by the 633 nm line of a He-Ne laser or the 647 nm of a Kr-Ar laser. Oyster[®]-650 is spectrally similar to Cy5[™] or Alexa Fluor[®] 647.

Cy3[™] and Cy5[™] are trademarks of GE Healthcare, Alexa Fluor[®] is a registered trademark of Invitrogen Inc.

Several of Denovo Biolabels products and product applications are covered by German and foreign patents and patents pending. No licensing agreement is necessary for use of Oyster[®]-500, -556, -645 and -656 as long as our products are utilized to produce another value added product for research or resale. For resale of our unmodified products as well as for all products involving Oyster[®]-550 or -650, a specific agreement or licensing agreement from Denovo Biolabels is necessary. All names containing the designation ® are registered with the German Patent and Trademark Office.

Copyright 2005, Denovo Biolabels GmbH. All rights reserved. This information is subject to change without notice.

For further information please contact our technical support

Denovo Biolabels GmbH

Mendelstr. 7

48149 Muenster, Germany

Phone +49 (0)251 980 2918

Fax +49 (0)251 980 2917

or visit our website: www.biolabels.com info@biolabels.com